

BenQ T3300

彩色超音波掃瞄器



128GB SSD
USB 3.0插槽
HDMI輸出及DICOM 3.0
1.5小時連續掃瞄
輕巧 堅固2.7kg
30秒快速開機
13.3" 多點觸控螢幕、易清潔、容易消毒



BenQ T3300

彩色超音波掃描器


是專門為定點照護(POC)所設計

藉由輕巧的機型和卓越的影像


讓您可以隨時隨地為您的病人提供醫療服務

臨床應用範圍


Abdominal, Small Parts, MSK, Vascular, Anesthesiology, Cardiology, OB/GYN.



Carotid Artery



Vascular measurement



Thyroid



Liver

C62B Convex Array

Frequency: 2~6 MHz

Application: Abdominal, MSK, OB/GYN, Fetal



LI154BH Linear Array

Frequency: 4~15 MHz

Application: Breast, Thyroid, Vascular, MSK, Nerve, Testicles



P42B6 Phased Array

Frequency: 2~4 MHz

Application: Cardiac, Abdominal



- Color Doppler / Power Doppler
- PW Doppler / CW Doppler
- Tissue Harmonic Imaging
- Compounding Imaging
- Q-Scan
- ENV (Enhanced Needle Visualization)

明基三豐醫療器材股份有限公司

台北總公司

台北市內湖區洲子街46號7樓

電話: +886-2-8797-5533

傳真: +886-2-8797-5188

www.BenQMedicalTech.com

info@BenQMedicalTech.com

Catalog No.: UL-VI707-tdr,
Copyright by BenQ Medical Technology Corp.

