



Make Veterinary Scanning Become Simple

## BenQ T3300v

- **Intuitive**
  - **Simplicity**
  - **High Mobility**
  - **Fast Boot-up**
  - **Easy to Clean**
- Multi Touch User Interface with Smart Finger Gestures
  - One-key Auto Optimization and Customized function keys
  - Portable Ultrasound can be Used Anywhere, Anytime
  - Less than 30 sec
  - Splash-Proof Screen with Virtual Keyboard



## Simple, Convenient, Quick to Diagnose


BenQ T3300v provides a friendly operating platform with smart finger gestures and multi touch user interface. Even if being a lightweight portable ultrasound, BenQ T3300v contains all necessary functions for a variety of situations. Relying on our outstanding image quality to aid fast and accurate diagnosis.

## Professional Design for Veterinary Use

- Animal ID info
- Specific Body Marks
- Figure Preset
- Measurement and Report

## Clinical Versatility with High Quality Image:

Never Compromise on Image Quality



### BenQ Medical Technology Corp.

Headquarter  
7FL 46 Zhou-z St., Neihs Dist.,  
Taipei 11493, Taiwan

Factory  
3F, No. 159, Shan-Ying Road, Gueishan District,  
Taoyuan City 333, Taiwan





Tel +886-2-8797-5533

Fax +886-2-8797-5188

[www.BenQMedicalTech.com](http://www.BenQMedicalTech.com)

[info@BenQMedicalTech.com](mailto:info@BenQMedicalTech.com)

Catalog No.: UL-V1909  
Copyright by BenQ Medical Technology Corp.

 4~15 MHz <b>Linear</b>	 2~6 MHz <b>Curved</b>	 3~8 MHz <b>Curved</b>	 3~8 MHz <b>Phased</b>
Application	Application	Application	Application
Small parts MSK Vein Kidney	Abdominal Pregnancy	Abdominal Pregnancy	Cardiology