BenQ H1300 Diagnostic Ultrasound System Image for Anywhere, Anytime



Mobile: Hand-held Ultrasound with 8" Multi-touch Screen

Light: 1.2 kg weight Quick: Boot-up in 30 sec

Consecutive: Rechargeable & swappable Li-ion battery

Reliable: 128 GB eMMC

Flexible: USB Type C, Mini-HDMI, DICOM 3.0, WiFi



Bring BenQ H1300 to the bedside

Can be used in various environments

- · ICU · ER · Ward · Field · Home care

Clinical versatility with High Quality Image:



Kidney

C62B_H1300 2~6 MHz

L154BH_H1300 P42B6_H1300



Vascular measurement



Curved

Application Abdomen OB/GYN Renal Urology MSK Chest

4~12 MHz Linear

Application Vascular Thyroid Breast MSK Nerve Chest

2~4 MHz **Phased**

Application Cardiac Chest

- Color Doppler / Power Doppler
- PW Doppler / CW Doppler
- Tissue Harmonic Imaging
- Compounding Imaging
- Q-Scan
- ENV (Enhanced Needle Visualization)

BenQ Medical Technology Corp.

Headquarter

7FL 46 Zhou-z St., Neihu Dist., Taipei 11493, Taiwan

Factory

3F, No. 159, Shan-Ying Road, Gueishan Hsiang, Taoyuan County 333, Taiwan, R.O.C.

Tel +886-2-8797-5533 Fax +886-2-8797-5188 www.BenQMedicalTech.com info@BenQMedicalTech.com